

P 20 dental dental plaster

| Properties | Special plaster of natural origin. (Plaster of Paris DIN EN ISO 6873 Typ 2) Can be coloured if necessary. | |
|----------------|--|---|
| Where to use | Basic material for dental laboratories. Hard and highly resistant universal dental impression compound. | |
| Advantages | Smooth surface Fast setting | |
| Technical data | Water/plaster ratio: Pot life (minutes): Setting time (minutes): Compressive strength (N/mm ²): Linear expansion (%): Hardness rating: Colour: | 0.58-0.60 7-9 11-13 > 50 < 0.15 Class 2 white |
| Packaging | Available in 25 kg paper bags, 40 bags/pallet, in heat-shrink wrapping in 5 kg paper bags, 200 bags/pallet, in heat-shrink wrapping | |
| Storage | In a covered place, in dry conditions, for approximately 6 months from the date of manufacture, in unopened bags Date of manufacture: see label on the bag. | |
| Manufacturer | CASEA GmbH., Germany | |
| Distributor | Széria-Trans Kereskedelmi és Szolgáltató Kft. H-1116 Budapest, Adony u. 7, Phone: +36-1-466-7758 H-5300 Karcag, Tilalmasi út 1, Phone: +36-59-311-333 | |

This data sheet was prepared to the best of our understanding. The data contained therein reflect our latest technical knowledge and observations. However, we shall not be held liable for the statutory provision of contents. Users shall honour any patent rights, and shall comply with the applicable laws and regulations at their own responsibility.